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- (1) Reddy ES, et al. The erg gene: a human gene related to the ets oncogene. Proc Natl Acad Sci U S A. 1987 Sep;84(17):6131-5.
- (2) Oikawa T. ETS transcription factors: possible targets for cancer therapy. Cancer Sci. 2004 Aug;95(8):626-33.
- (3) FitzGerald LM, et al. Association of TMPRSS2-ERG gene fusion with clinical characteristics and outcomes: results from a population-based study of prostate cancer. BMC Cancer. 2008 Aug 11;8:230.
- (4) Kumar-Sinha C, et al. Recurrent gene fusions in prostate cancer. Nat Rev Cancer. 2008 Jul;8(7):497-511.
- (5) Leong M, et al. Overexpression of truncated ERG from TMPRSS2-ERG fusion and prostate cancer development. Pathol Lab Med Int 2009;1:13-21
- (6) Furusato B, et al. ERG oncoprotein expression in prostate cancer: clonal progression of ERG-positive tumor cells and potential for ERG-based stratification. Prostate Cancer Prostatic Dis. 2010 Sep;13(3):228-37.
- (7) Minner S, et al. ERG status is unrelated to PSA recurrence in radically operated prostate cancer in the absence of antihormonal therapy. Clin Cancer Res. 2011 Sep 15;17(18):5878-88.
- (8) Falzarano SM, et al. ERG gene rearrangement status in prostate cancer detected by immunohistochemistry. Virchows Arch. 2011 Oct;459(4):441-7.
- (9) Chaux A, et al. Immunohistochemistry for ERG expression as a surrogate for TMPRSS2-ERG fusion detection in prostatic adenocarcinomas. Am J Surg Pathol. 2011 Jul;35(7):1014-20.
- (10) Park K, et al. Antibody-based detection of ERG rearrangement-positive prostate cancer. Neoplasia. 2010 Jul;12(7):590-8.
- (11) Huang L, et al. The oncogenic gene fusion TMPRSS2: ERG is not a diagnostic or prognostic marker for ovarian cancer. Int J Clin Exp Pathol. 2011;4(7):644-50. Epub 2011 Sep 15.
- (12) Miettinen M, et al. ERG transcription factor as an immunohistochemical marker for vascular endothelial tumors and prostatic carcinoma. Am J Surg Pathol. 2011 Mar;35(3):432-41.
- (13) Rosen P, et al. Clinical potential of the ERG oncoprotein in prostate cancer. Nat Rev Urol. 2012 Feb 14;9(3):131-7.
- (14) Braun M, et al. ERG protein expression and genomic rearrangement status in primary and metastatic prostate cancer—a comparative study of two monoclonal antibodies. Prostate Cancer Prostatic Dis. 2012 Jun;15(2):165-9.
- (15) Hoogland AM, et al. ERG immunohistochemistry is not predictive for PSA recurrence, local recurrence or overall survival after radical prostatectomy for prostate cancer. Mod Pathol. 2012 Mar;25(3):471-9.

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