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- (1) Cunningham BA, et al. Neural cell adhesion molecule: structure, immunoglobulin-like domains, cell surface modulation, and alternative RNA splicing. Science. 1987 May 15;236(4803):799-806.
- (2) Mechtersheimer G, et al. Expression of the natural killer cell-associated antigens CD56 and CD57 in human neural and striated muscle cells and in their tumors. Cancer Res. 1991 Feb 15;51(4):1300-7.
- (3) Kibbelaar RE, et al. Neural cell adhesion molecule expression, neuroendocrine differentiation and prognosis in lung carcinoma. Eur J Cancer. 1991;27(4):431-5.
- (4) Michalides R, et al. NCAM and lung cancer. Int J Cancer Suppl. 1994;8:34-7.
- (5) Gerardy-Schahn R, et al. Hot spots of antigenicity in the neural cell adhesion molecule NCAM. Int J Cancer Suppl. 1994;8:38-42.
- (6) Kaufmann O, et al. M. Utility of 123C3 monoclonal antibody against CD56 (NCAM) for the diagnosis of small cell carcinomas on paraffin sections. Hum Pathol. 1997 Dec;28(12):1373-8.
- (7) Trejo O, et al. Atypical cells in human cutaneous re-excision scars for melanoma express p75NGFR, C56/N-CAM and GAP-43: evidence of early Schwann cell differentiation. J Cutan Pathol. 2002 Aug;29(7):397-406.
- (8) Ely SA, et al. Expression of CD56/neural cell adhesion molecule correlates with the presence of lytic bone lesions in multiple myeloma and distinguishes myeloma from monoclonal gammopathy of undetermined significance and lymphomas with plasmacytoid differentiation. Am J Pathol. 2002 Apr;160(4):1293-9.
- (9) Tao J, et al. Aggressive Epstein-Barr virus-associated, CD8+, CD30+, CD56+, surface CD3-, natural killer (NK)-like cytotoxic T-cell lymphoma. Am J Surg Pathol. 2002 Jan;26(1):111-8.
- (10) Sumi M, et al. Natural killer cell lymphoma in the duodenum. Leuk Lymphoma. 2003 Jan;44(1):201-4.

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