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- (1) Liu H, et al. Immunohistochemical evaluation of GATA3 expression in tumors and normal tissues: a useful immunomarker for breast and urothelial carcinomas. Am J Clin Pathol. 2012 Jul;138(1):57-64.
- (2) Chang A, et al. Utility of GATA3 immunohistochemistry in differentiating urothelial carcinoma from prostate adenocarcinoma and squamous cell carcinomas of the uterine cervix, anus, and lung. Am J Surg Pathol. 2012 Oct;36(10):1472-6.
- (3) Ordóñez NG. Value of GATA3 immunostaining in tumor diagnosis: a review. Adv Anat Pathol. 2013 Sep;20(5):352-60.
- (4) Gonzalez-Roibon N, et al. The role of GATA binding protein 3 in the differential diagnosis of collecting duct and upper tract urothelial carcinomas. Hum Pathol. 2013 Dec;44(12):2651-7.
- (5) Miettinen M, et al. GATA3: a multispecific but potentially useful marker in surgical pathology: a systematic analysis of 2500 epithelial and nonepithelial tumors. Am J Surg Pathol. 2014 Jan;38(1):13-22.
- (6) Gonzalez-Roibon N, et al. Comprehensive profile of GATA binding protein 3 immunohistochemical expression in primary and metastatic renal neoplasms. Hum Pathol. 2014 Feb;45(2):244-8.
- (7) Shield PW, et al. GATA3: a promising marker for metastatic breast carcinoma in serous effusion specimens. Cancer Cytopathol. 2014 Apr;122(4):307-12.
- (8) Liu H, et al. Immunohistochemical evaluation of GATA-3 expression in ER-negative breast carcinomas. Am J Clin Pathol. 2014 May;141(5):648-55.
- (9) Clark BZ, et al. GATA-3 immunoreactivity in breast, bladder, gynecologic tract, and other cytokeratin 7-positive carcinomas. Am J Clin Pathol. 2014 Jul;142(1):64-71.
- (10) Liang Y, et al. Differential expression of GATA-3 in urothelial carcinoma variants. Hum Pathol. 2014 Jul;45(7):1466-72.
- (11) Krings G, et al. Diagnostic utility and sensitivities of GATA3 antibodies in triple-negative breast cancer. Hum Pathol. 2014 Nov;45(11):2225-32.
- (12) Banet N, et al. GATA-3 expression in trophoblastic tissues: an immunohistochemical study of 445 cases, including diagnostic utility. Am J Surg Pathol. 2015 Jan;39(1):101-8.
- (13) Sangoi AR, et al. The Novel Marker GATA3 is Significantly More Sensitive Than Traditional Markers Mammaglobin and GCDFP15 for Identifying Breast Cancer in Surgical and Cytology Specimens of Metastatic and Matched Primary Tumors. Appl Immunohistochem Mol Morphol. 2016 Apr;24(4):229-37.

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