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- (1) Kimura N, et al. Immunocytochemical localization of androgen receptor with polyclonal antibody in paraffin-embedded human tissues. J Histochem Cytochem. 1993 May;41(5):671-8.
- (2) Chang C, et al. Androgen receptor: an overview. Crit Rev Eukaryot Gene Expr. 1995;5(2):97-125.
- (3) Ozanne DM, et al. Androgen receptor nuclear translocation is facilitated by the f-actin cross-linking protein filamin. Mol Endocrinol. 2000 Oct;14(10):1618-26.
- (4) Niemeier LA, et al. Androgen receptor in breast cancer: expression in estrogen receptor-positive tumors and in estrogen receptor-negative tumors with apocrine differentiation. Mod Pathol. 2010 Feb;23(2):205-12.
- (5) Nodin B, et al. Increased androgen receptor expression in serous carcinoma of the ovary is associated with an improved survival. J Ovarian Res. 2010 Jun 17;3:14.
- (6) Collins LC, et al. Androgen receptor expression in breast cancer in relation to molecular phenotype: results from the Nurses' Health Study. Mod Pathol. 2011 Jul;24(7):924-31.
- (7) Yu Q, et al. Expression of androgen receptor in breast cancer and its significance as a prognostic factor. Ann Oncol. 2011 Jun;22(6):1288-1294.
- (8) Tsutsumi Y. Apocrine carcinoma as triple-negative breast cancer: novel definition of apocrine-type carcinoma as estrogen/progesterone receptor-negative and androgen receptor-positive invasive ductal carcinoma. Jpn J Clin Oncol. 2012 May;42(5):375-86.
- (9) Gucalp A, et al. Translational Breast Cancer Research Consortium (TBCRC 011). Phase II trial of bicalutamide in patients with androgen receptor-positive, estrogen receptor-negative metastatic Breast Cancer. Clin Cancer Res. 2013 Oct 1;19(19):5505-12.
- (10) Safarpour D, et al. Androgen receptor (AR) expression in 400 breast carcinomas: is routine AR assessment justified? Am J Cancer Res. 2014 Jul 16;4(4):353-68.
- (11) Queisser A, et al. Comparison of different prostatic markers in lymph node and distant metastases of prostate cancer. Mod Pathol. 2015 Jan;28(1):138-45.
- (12) Schalken J, et al. Enzalutamide: targeting the androgen signalling pathway in metastatic castration-resistant prostate cancer. BJU Int. 2016 Feb;117(2):215-25.

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