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- (1) Nielsen GP, et al. Immunohistochemical survey of p16INK4A expression in normal human adult and infant tissues. Lab Invest. 1999 Sep;79(9):1137-43.
- (2) Romagosa C, et al. p16(Ink4a) overexpression in cancer: a tumor suppressor gene associated with senescence and high-grade tumors. Oncogene. 2011 May 5;30(18):2087-97.
- (3) Jedpiyawongse A, et al. Immunohistochemical overexpression of p16 protein associated with cervical cancer in Thailand. Asian Pac J Cancer Prev. 2008 Oct-Dec;9(4):625-30.
- (4) Shelton J, et al. p16 immunohistochemistry in oropharyngeal squamous cell carcinoma: a comparison of antibody clones using patient outcomes and high-risk human papillomavirus RNA status. Mod Pathol. 2017 Sep;30(9):1194-1203.
- (5) Yang CH, et al. Expressions of p16 and p27 in urothelial carcinoma and their prognostic value. Kaohsiung J Med Sci. 2014 Sep;30(9):453-8.
- (6) Al-Ahwal M, et al. p16 protein is upregulated in a stepwise fashion in colorectal adenoma and colorectal carcinoma. Saudi J Gastroenterol. 2016 Nov;22(6):435-440.
- (7) Stoler MH, et al. Routine Use of Adjunctive p16 Immunohistochemistry Improves Diagnostic Agreement of Cervical Biopsy Interpretation: Results From the CERTAIN Study. Am J Surg Pathol. 2018 Aug;42(8):1001-1009.
- (8) Darragh TM, et al. The Lower Anogenital Squamous Terminology Standardization Project for HPV-Associated Lesions: background and consensus recommendations from the College of American Pathologists and the American Society for Colposcopy and Cervical Pathology. Arch Pathol Lab Med. 2012 Oct;136(10):1266-97.
- (9) Kalof AN, et al. Our approach to squamous intraepithelial lesions of the uterine cervix. J Clin Pathol. 2007 May;60(5):449-55.

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