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- (1) Bafna S, et al. Membrane-bound mucins: the mechanistic basis for alterations in the growth and survival of cancer cells. *Oncogene*. 2010 May 20;29(20):2893-904.
- (2) Lau SK, et al. Differential expression of MUC1, MUC2, and MUC5AC in carcinomas of various sites: an immunohistochemical study. *Am J Clin Pathol*. 2004 Jul;122(1):61-9.
- (3) Walsh MD, et al. Expression of MUC2, MUC5AC, MUC5B, and MUC6 mucins in colorectal cancers and their association with the CpG island methylator phenotype. *Mod Pathol*. 2013 Dec;26(12):1642-56.
- (4) Sasaki M, et al. Expression Profiles of MUC Mucin Core Protein in the Intrahepatic Biliary System: Physiological Distribution and Pathological Significance. *Acta Histochem Cytochem*. 2005;38(5):295-303.
- (5) Nakajima K, et al. Expression of gastric gland mucous cell-type mucin in normal and neoplastic human tissues. *J Histochem Cytochem*. 2003 Dec;51(12):1689-98.
- (6) Arul GS, et al. Mucin gene expression in Barrett's oesophagus: an in situ hybridisation and immunohistochemical study. *Gut*. 2000 Dec;47(6):753-61.
- (7) Tsukashita S, et al. MUC gene expression and histogenesis of adenocarcinoma of the stomach. *Int J Cancer*. 2001 Oct 15;94(2):166-70.
- (8) Lee HS, et al. MUC1, MUC2, MUC5AC, and MUC6 expressions in gastric carcinomas: their roles as prognostic indicators. *Cancer*. 2001 Sep 15;92(6):1427-34.
- (9) Betge J, et al. MUC1, MUC2, MUC5AC, and MUC6 in colorectal cancer: expression profiles and clinical significance. *Virchows Arch*. 2016 Sep;469(3):255-65.
- (10) Kobayashi M, et al. Reappraisal of the Immunophenotype of Pancreatic Intraductal Papillary Mucinous Neoplasms (IPMNs)-Gastric Pyloric and Small Intestinal Immunophenotype Expression in Gastric and Intestinal Type IPMNs-. *Acta Histochem Cytochem*. 2014 May 1;47(2):45-57.

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