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- (1)Kellner S, et al. Transcriptional regulation of the Oct4 gene, a master gene for pluripotency. *Histol Histopathol*. 2010 Mar;25(3):405-12.
- (2)Rizzino A. Sox2 and Oct-3/4: a versatile pair of master regulators that orchestrate the self-renewal and pluripotency of embryonic stem cells. *Wiley Interdiscip Rev Syst Biol Med*. 2009 Sep-Oct;1(2):228-236.
- (3)Hattab EM, Tu PH, Wilson JD, Cheng L. OCT4 immunohistochemistry is superior to placental alkaline phosphatase (PLAP) in the diagnosis of central nervous system germinoma. *Am J Surg Pathol*. 2005 Mar;29(3):368-71.
- (4)Looijenga LH, et al. POU5F1 (OCT3/4) identifies cells with pluripotent potential in human germ cell tumors. *Cancer Res*. 2003 May 1;63(9):2244-50.
- (5)Rajpert-De Meyts E, et al. Diagnostic markers for germ cell neoplasms: from placental-like alkaline phosphatase to micro-RNAs. *Folia Histochem Cytobiol*. 2015;53(3):177-88.
- (6)Rajpert-De Meyts E, et al. Developmental expression of POU5F1 (OCT3/4) in normal and dysgenetic human gonads. *Hum Reprod*. 2004 Jun;19(6):1338-44.
- (7)Ulbricht TM, et al: Members of the ISUP Immunohistochemistry in Diagnostic Urologic Pathology Group. Best practices recommendations in the application of immunohistochemistry in testicular tumors: report from the International Society of Urological Pathology consensus conference. *Am J Surg Pathol*. 2014 Aug;38(8):e50-9.
- (8)Looijenga LHJ, et al. Report From the International Society of Urological Pathology (ISUP) Consultation Conference on Molecular Pathology of Urogenital Cancers: IV: Current and Future Utilization of Molecular-Genetic Tests for Testicular Germ Cell Tumors. *Am J Surg Pathol*. 2020 Jul;44(7):e66-e79.
- (9)Già N, et al. Ovarian Germ Cell Tumors: Pictorial Essay. *Diagnostics (Basel)*. 2022 Aug 24;12(9):2050.
- (10)Mei K, et al. Diagnostic utility of SALL4 in primary germ cell tumors of the central nervous system: a study of 77 cases. *Mod Pathol*. 2009 Dec;22(12):1628-36.

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