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- (1) Matsuda R, et al. Expression of the c-kit protein in human solid tumors and in corresponding fetal and adult normal tissues. *Am J Pathol.* 1993 Jan;142(1):339-46.
- (2) Tsuura Y, et al. Preferential localization of c-kit product in tissue mast cells, basal cells of skin, epithelial cells of breast, small cell lung carcinoma and seminoma/dysgerminoma in human: immunohistochemical study on formalin-fixed, paraffin-embedded tissues. *Virchows Arch.* 1994;424(2):135-41.
- (3) Yamataka A, et al. A lack of intestinal pacemaker (c-kit) in aganglionic bowel of patients with Hirschsprung's disease. *J Pediatr Surg.* 1995 Mar;30(3):441-4.
- (4) Yamataka A, et al. Lack of intestinal pacemaker (C-KIT-positive) cells in infantile hypertrophic pyloric stenosis. *J Pediatr Surg.* 1996 Jan;31(1):96-8; discussion 98-9.
- (5) Yamataka A, et al. Localization of intestinal pacemaker cells and synapses in the muscle layers of a patient with colonic hypoganglionosis. *J Pediatr Surg.* 1996 Apr;31(4):584-7.
- (6) Yamataka A, et al. Intestinal pacemaker C-KIT+ cells and synapses in allied Hirschsprung's disorders. *J Pediatr Surg.* 1997 Jul;32(7):1069-74.
- (7) Arber DA, et al. Weiss LM. Paraffin section detection of the c-kit gene product (CD117) in human tissues: value in the diagnosis of mast cell disorders. *Hum Pathol.* 1998 May;29(5):498-504.
- (8) Kindblom LG, et al. Gastrointestinal pacemaker cell tumor (GIPACT): gastrointestinal stromal tumors show phenotypic characteristics of the interstitial cells of Cajal. *Am J Pathol.* 1998 May;152(5):1259-69.
- (9) Yamataka A, et al. Abnormal distribution of intestinal pacemaker (C-KIT-positive) cells in an infant with chronic idiopathic intestinal pseudoobstruction. *J Pediatr Surg.* 1998 Jun;33(6):859-62.
- (10) Hirota S, et al. Gain-of-function mutations of c-kit in human gastrointestinal stromal tumors. *Science.* 1998 Jan 23;279(5350):577-80.
- (11) Ernst SI, et al. KIT mutation portends poor prognosis in gastrointestinal stromal/smooth muscle tumors. *Lab Invest.* 1998 Dec;78(12):1633-6.
- (12) Komuro T, et al. Ultrastructural characterization of the interstitial cells of Cajal. *Arch Histol Cytol.* 1999 Oct;62(4):295-316.
- (13) 秋濱 玄 他. Gastrointestinal stromal tumor の組織発生及び悪性度の解析. *岩手医誌.* 1999;51:381-390.
- (14) Natkunam Y, et al. Utility of paraffin section immunohistochemistry for C-KIT (CD117) in the differential diagnosis of systemic mast cell disease involving the bone marrow. *Am J Surg Pathol.* 2000 Jan;24(1):81-91.
- (15) Maeyama H, et al. Familial gastrointestinal stromal tumor with hyperpigmentation: association with a germline mutation of the c-kit gene. *Gastroenterology.* 2001 Jan;120(1):210-5.

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