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- (1) Lo Coco F, et al. Rearrangements of the BCL6 gene in diffuse large cell non-Hodgkin's lymphoma. Blood. 1994 Apr 1;83(7):1757-9.
- (2) Pescarmona E, et al. Analysis of the BCL-6 gene configuration in diffuse B-cell non-Hodgkin's lymphomas and Hodgkin's disease. J Pathol. 1995 Sep;177(1):21-5.
- (3) Otsuki T, et al. Analysis of LAZ3 (BCL-6) status in B-cell non-Hodgkin's lymphomas: results of rearrangement and gene expression studies and a mutational analysis of coding region sequences. Blood. 1995 May 15;85(10):2877-84.
- (4) Flenghi L, et al. Monoclonal antibodies PG-B6a and PG-B6p recognize, respectively, a highly conserved and a formol-resistant epitope on the human BCL-6 protein amino-terminal region. Am J Pathol. 1996 May;148(5):1543-55.
- (5) Yang B, et al. Clinicopathologic reassessment of primary cutaneous B-cell lymphomas with immunophenotypic and molecular genetic characterization. Am J Surg Pathol. 2000 May;24(5):694-702.
- (6) Dunphy CH, et al. Paraffin immunoreactivity of CD10, CDw75, and Bcl-6 in follicle center cell lymphoma. Leuk Lymphoma. 2001 May;41(5-6):585-92.
- (7) Fernández-Vázquez A, et al. Primary cutaneous large B-cell lymphoma: the relation between morphology, clinical presentation, immunohistochemical markers, and survival. Am J Surg Pathol. 2001 Mar;25(3):307-15.
- (8) Hans CP, et al. Confirmation of the molecular classification of diffuse large B-cell lymphoma by immunohistochemistry using a tissue microarray. Blood. 2004 Jan 1;103(1):275-82.

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