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- (1) Lan MS, et al. Genomic organization, 5'-upstream sequence, and chromosomal localization of an insulinoma-associated intronless gene, IA-1. J Biol Chem. 1994 May 13;269(19):14170-4.
- (2) Goto Y, et al. A novel human insulinoma-associated cDNA, IA-1, encodes a protein with "zinc-finger" DNA-binding motifs. J Biol Chem. 1992 Jul 25;267(21):15252-7.
- (3) Chen C, et al. Insulinoma-Associated-1: From Neuroendocrine Tumor Marker to Cancer Therapeutics. Mol Cancer Res. 2019 Aug;17(8):1597-1604.
- (4) Gierl MS, et al. The zinc-finger factor Insm1 (IA-1) is essential for the development of pancreatic beta cells and intestinal endocrine cells. Genes Dev. 2006 Sep 1;20(17):2465-78.
- (5) Yoshida A, et al. INSM1 expression and its diagnostic significance in extraskeletal myxoid chondrosarcoma. Mod Pathol. 2018 May;31(5):744-752.
- (6) Rosenbaum JN, et al. INSM1: A Novel Immunohistochemical and Molecular Marker for Neuroendocrine and Neuroepithelial Neoplasms. Am J Clin Pathol. 2015 Oct;144(4):579-91.
- (7) McHugh KE, et al. INSM1 Is a Highly Specific Marker of Neuroendocrine Differentiation in Primary Neoplasms of the Gastrointestinal Tract, Appendix, and Pancreas. Am J Clin Pathol. 2020 May 5;153(6):811-820.
- (8) Mukhopadhyay S, et al. Insulinoma-associated protein 1 (INSM1) is a sensitive and highly specific marker of neuroendocrine differentiation in primary lung neoplasms: an immunohistochemical study of 345 cases, including 292 whole-tissue sections. Mod Pathol. 2019 Jan;32(1):100-109.
- (9) Kuji S, et al. Immunosenitivity and specificity of insulinoma-associated protein 1 (INSM1) for neuroendocrine neoplasms of the uterine cervix. J Gynecol Oncol. 2023 Jan;34(1):e1.
- (10) Metovic J, et al. INSM1 Expression in Breast Neoplasms with Neuroendocrine Features. Endocr Pathol. 2021 Dec;32(4):452-460.
- (11) Rooper LM, et al. INSM1 Demonstrates Superior Performance to the Individual and Combined Use of Synaptophysin, Chromogranin and CD56 for Diagnosing Neuroendocrine Tumors of the Thoracic Cavity. Am J Surg Pathol. 2017 Nov;41(11):1561-1569.

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