NICHIREI BIOSCIENCES INC.

Focus Report vol.3

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Contents

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
- 1) Hillerdal, G.: Malignant mesothelioma 1982: review of 4710 published cases. *Br J Dis Chest* 1983, 77: 321-343
- 2) Neumann, V., Gunthe, S., Mülle, K. M., et al.: Malignant mesothelioma--German mesothelioma register 1987-1999. *Int Arch Occup Environ Health* 2001, 74: 383-395
- 3) Flores, R. M., Pass, H. I., Seshan, V. E., et al.: Extrapleural pneumonectomy versus pleurectomy/decortication in the surgical management of malignant pleural mesothelioma: results in 663 patients. *J Thorac Cardiovasc Surg* 2008, 135: 620-626, 626 e621-623
- 4) Ettinger, D. S., Wood, D. E., Akerley, W., et al. (2015). Malignant pleural mesothelioma *NCCN Clinical Practice Guidelines in Oncology* National Comprehensive Cancer Network, Inc.
- 5) Metintas, M., Ak, G., Dundar, E., et al.: Medical thoracoscopy vs CT scan-guided Abrams pleural needle biopsy for diagnosis of patients with pleural effusions: a randomized, controlled trial. *Chest* 2010, 137: 1362-1368
- 6) Attanoos, R. L., & Gibbs, A. R.: 'Pseudomesotheliomatous' carcinomas of the pleura: a 10-year analysis of cases from the Environmental Lung Disease Research Group, Cardiff. *Histopathology* 2003, 43: 444-452
- 7) Harwood, T. R., Gracey, D. R., & Yokoo, H.: Pseudomesotheliomatous carcinoma of the lung. A variant of peripheral lung cancer. *Am J Clin Pathol* 1976, 65: 159-167

- 8) Husain, A. N., Colby, T. V., Ordonez, N. G., et al.: Guidelines for Pathologic Diagnosis of Malignant Mesothelioma 2017 Update of the Consensus Statement From the International Mesothelioma Interest Group. *Archives of Pathology and Laboratory Medicine* 2018, 142: 89-108
- 9) Ordonez, N. G.: Value of claudin-4 immunostaining in the diagnosis of mesothelioma. *Am J Clin Pathol* 2013, 139: 611-619
- 10) Ordonez, N. G.: The diagnostic utility of immunohistochemistry in distinguishing between mesothelioma and renal cell carcinoma: a comparative study. *Hum Pathol* 2004, 35: 697-710
- 11) Butnor, K. J., Nicholson, A. G., Allred, D. C., et al.: Expression of renal cell carcinoma-associated markers erythropoietin, CD10, and renal cell carcinoma marker in diffuse malignant mesothelioma and metastatic renal cell carcinoma. *Arch Pathol Lab Med* 2006, 130: 823-827
- 12) Miettinen, M., McCue, P. A., Sarlomo-Rikala, M., et al.: GATA3: a multispecific but potentially useful marker in surgical pathology: a systematic analysis of 2500 epithelial and nonepithelial tumors. *Am J Surg Pathol* 2014, 38: 13-22
- 13) Chapel, D. B., Husain, A. N., Krausz, T., et al.: PAX8 Expression in a Subset of Malignant Peritoneal Mesotheliomas and Benign Mesothelium has Diagnostic Implications in the Differential Diagnosis of Ovarian Serous Carcinoma. *Am J Surg Pathol* 2017, 41: 1675-1682
- 14) Wu, D., Hiroshima, K., Yusa, T., et al.: Usefulness of p16/CDKN2A fluorescence in situ hybridization and BAP1 immunohistochemistry for the diagnosis of biphasic mesothelioma. *Ann Diagn Pathol* 2017, 26: 31-37
- 15) Wu, D., Hiroshima, K., Matsumoto, S., et al.: Diagnostic usefulness of p16/CDKN2A FISH in distinguishing between sarcomatoid mesothelioma and fibrous pleuritis. *Am J Clin Pathol* 2013, 139: 39-46
- 16) 環境省：平成28年度石綿関連疾患に係る医学的所見の解析調査業務（FISH法を用いた中皮腫診断法の開発に関する調査編）報告書（主任研究者：廣島健三）。（2016）。
- 17) Tochigi, N., Attanoos, R., Chirieac, L. R., et al.: p16 Deletion in sarcomatoid tumors of the lung and pleura. *Arch Pathol Lab Med* 2013, 137: 632-636
- 18) Hwang, H. C., Pyott, S., Rodriguez, S., et al.: BAP1 Immunohistochemistry and p16 FISH in the Diagnosis of Sarcomatous and Desmoplastic Mesotheliomas. *Am J Surg Pathol* 2016, 40: 714-718
- 19) Husain, A. N., Colby, T., Ordonez, N., et al.: Guidelines for pathologic diagnosis of malignant mesothelioma: 2012 update of the consensus statement from the International Mesothelioma Interest Group. *Arch Pathol Lab Med* 2013, 137: 647-667
- 20) Churg, A., Cagle, P. T., & Roggli, V. L. (2006). *Tumors of the Serosal Membrane*. Washington, DC: ARP Press.
- 21) Berg, K. B., & Churg, A.: GATA3 Immunohistochemistry for Distinguishing Sarcomatoid and Desmoplastic Mesothelioma From Sarcomatoid Carcinoma of the Lung. *Am J Surg Pathol* 2017, 41: 1221-1225
- 22) Amatya, V. J., Kushitani, K., Mawas, A. S., et al.: MUC4, a novel immunohistochemical marker identified by gene expression profiling, differentiates pleural sarcomatoid mesothelioma from lung sarcomatoid carcinoma. *Mod Pathol* 2017, 30: 672-681
- 23) Berg, K. B., Ionescu, D., & Churg, A. M.: MUC4 Staining in Sarcomatoid Carcinomas. *Mod Pathol* 2019, 32: 157
- 24) Terra, S., Roden, A., Yi, E., et al.: Utility of immunohistochemistry for MUC4 and GATA3 to aid in the distinction of sarcomatoid mesothelioma from sarcomatoid carcinoma. *Mod Pathol Supplement* 2. 2019, 32: 59
- 25) Tsuji, S., Washimi, K., Kageyama, T., et al.: HEG1 is a novel mucin-like membrane protein that serves as a diagnostic and therapeutic target for malignant mesothelioma. *Sci Rep* 2017, 7: 45768
- 26) Hjerpe, A., Ascoli, V., Bedrossian, C. W., et al.: Guidelines for the cytopathologic diagnosis of epithelioid and mixed-type malignant mesothelioma. Complementary statement from the international mesothelioma interest group, also endorsed by the international academy of cytology and the papanicolaou society of cytopathology. *Acta Cytol* 2015, 59: 2-16
- 27) Hwang, H. C., Sheffield, B. S., Rodriguez, S., et al.: Utility of BAP1 immunohistochemistry and p16 (CDKN2A) FISH in the diagnosis of malignant mesothelioma in effusion cytology specimens. *Am J Surg Pathol* 2016, 40: 120-126
- 28) Walts, A. E., Hiroshima, K., McGregor, S. M., et al.: BAP1 immunostain and CDKN2A (p16) FISH analysis: clinical applicability for the diagnosis of malignant mesothelioma in effusions. *Diagn Cytopathol* 2016, 44: 599-606
- 29) Hiroshima, K., Wu, D., Hasegawa, M., et al.: Cytologic differential diagnosis of malignant mesothelioma and reactive mesothelial cells with FISH analysis of p16. *Diagn Cytopathol* 2016, 44: 591-598
- 30) Hasegawa, S., Kondo, N., Matsumoto, S., et al.: Practical approaches to diagnose and treat for T0 malignant pleural mesothelioma: a proposal for diagnostic total parietal pleurectomy. *Int J Clin Oncol* 2012, 17: 33-39
- 31) Whitaker, D., Henderson, D. W., & Shilkin, K. B.: The concept of mesothelioma in situ: implications for diagnosis and histogenesis. *Semin Diagn Pathol* 1992, 9: 151-161
- 32) Henderson, D. W., Shilkin, K. B., & Whitaker, D.: Reactive mesothelial hyperplasia vs mesothelioma, including mesothelioma in situ: a brief review. *Am J Clin Pathol* 1998, 110: 397-404
- 33) Churg, A., Hwang, H., Tan, L., et al.: Malignant mesothelioma in situ. *Histopathology* 2018, 72: 1033-1038
- 34) Scherpereel, A., Astoul, P., Baas, P., et al.: Guidelines of the European Respiratory Society and the European Society of Thoracic Surgeons for the management of malignant pleural mesothelioma. *Eur Respir J* 2010, 35: 479-495
- 35) BTS statement on malignant mesothelioma in the UK, 2007. *Thorax* 2007, 62 Suppl 2: ii1-ii19

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p16 HD, homozygous deletion of *p16* gene; FISH, fluorescence *in situ* hybridization.

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
- 1) Husain AN et al., Guidelines for pathologic diagnosis of malignant mesothelioma – 2017 update of the consensus statement from the International Mesothelioma Interest Group. Arch Pathol Lab Med 2018;142:89-108
- 2) Hida T et al., Immunohistochemical detection of MTAP and BAP1 protein loss for mesothelioma diagnosis: Comparison with 9p21 FISH and BAP1 immunohistochemistry. Lung Cancer. 2017; 104:98-105
- 3) Kinoshita Y et al., A combination of MTAP and BAP1 immunohistochemistry in pleural effusion cytology for the diagnosis of mesothelioma. Cancer Cytopathol. 2018;126:54-63

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