**NICHIREI BIOSCIENCES INC.**

**Focus Report vol.4**

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## Contents

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- 1) Rusch, V. W.: A proposed new international TNM staging system for malignant pleural mesothelioma. From the International Mesothelioma Interest Group. *Chest* 1995, 108: 1122-1128
- 2) Nowak, A. K., Chansky, K., Rice, D. C., et al.: The IASLC Mesothelioma Staging Project: Proposals for Revisions of the T Descriptors in the Forthcoming Eighth Edition of the TNM Classification for Pleural Mesothelioma. *J Thorac Oncol* 2016, 11: 2089-2099
- 3) Gill, R. R., Nowak, A. K., Giroux, D. J., et al.: The International Association for the Study of Lung Cancer Mesothelioma Staging Project: Proposals for Revisions of the "T" Descriptors in the Forthcoming Ninth Edition of the TNM Classification for Pleural Mesothelioma. *J Thorac Oncol* 2024, 19: 1310-1325
- 4) Hillerdal, G.: Malignant mesothelioma 1982: review of 4710 published cases. *Br J Dis Chest* 1983, 77: 321-343
- 5) Neumann, V., Gunthe, S., Mülle, K. M., et al.: Malignant mesothelioma--German mesothelioma register 1987-1999. *Int Arch Occup Environ Health* 2001, 74: 383-395
- 6) Flores, R. M., Pass, H. I., Seshan, V. E., et al.: Extrapleural pneumonectomy versus pleurectomy/decortication in the surgical management of malignant pleural mesothelioma: results in 663 patients. *J Thorac Cardiovasc Surg* 2008, 135: 620-626, 626 e621-623
- 7) Baas, P., Scherpereel, A., Nowak, A. K., et al.: First-line nivolumab plus ipilimumab in unresectable malignant pleural mesothelioma (CheckMate 743): a multicentre, randomised, open-label, phase 3 trial. *Lancet* 2021, 397: 375-386
- 8) Lim, E., Waller, D., Lau, K., et al.: Extended pleurectomy decortication and chemotherapy versus chemotherapy alone for pleural mesothelioma (MARS 2): a phase 3 randomised controlled trial. *Lancet Respir Med* 2024, 12: 457-466
- 9) Mesothelioma: Pleural, Version 2.2026 — October 3, 2025 NCCN Clinical Practice Guidelines in Oncology
- 10) Metintas, M., Ak, G., Dundar, E., et al.: Medical thoracoscopy vs CT scan-guided Abrams pleural needle biopsy for diagnosis of patients with pleural effusions: a randomized, controlled trial. *Chest* 2010, 137: 1362-1368
- 11) Attanoos, R. L., & Gibbs, A. R.: 'Pseudomesotheliomatous' carcinomas of the pleura: a 10-year analysis of cases from the Environmental Lung Disease Research Group, Cardiff. *Histopathology* 2003, 43: 444-452
- 12) Harwood, T. R., Gracey, D. R., & Yokoo, H.: Pseudomesotheliomatous carcinoma of the lung. A variant of peripheral lung cancer. *Am J Clin Pathol* 1976, 65: 159-167
- 13) Husain, A. N., Chapel, D. B., Attanoos, R., et al.: Guidelines for Pathologic Diagnosis of Mesothelioma: 2023 Update of the Consensus Statement From the International Mesothelioma Interest Group. *Arch Pathol Lab Med* 2024, 148: 1251-1271
- 14) Matsuda, M., Ninomiya, H., Wakejima, R., et al.: Calretinin-expressing lung adenocarcinoma: Distinct characteristics of advanced stages, smoker-type features, and rare expression of other mesothelial markers are useful to differentiate epithelioid mesothelioma. *Pathol Res Pract* 2020, 216: 152817
- 15) Tsuji, S., Washimi, K., Kageyama, T., et al.: HEG1 is a novel mucin-like membrane protein that serves as a diagnostic and therapeutic target for malignant mesothelioma. *Sci Rep* 2017, 7: 45768
- 16) Hiroshima K., Wu D., Koh E., et al.: Membranous HEG1 expression is a useful marker in the differential diagnosis of epithelioid and biphasic malignant mesothelioma versus carcinomas. *Pathol Int*. 2021, 71(9):604-613.
- 17) Hiroshima, K., Wu, D., Hamakawa, S., et al.: HEG1, BAP1, and MTAP are useful in cytologic diagnosis of malignant mesothelioma with effusion. *Diagn Cytopathol* 2021, 49: 622-632
- 18) Ordonez, N. G.: Value of claudin-4 immunostaining in the diagnosis of mesothelioma. *Am J Clin Pathol* 2013, 139: 611-619
- 19) Lonardi, S., Manera, C., Marucci, R., et al.: Usefulness of Claudin 4 in the cytological diagnosis of serosal effusions. *Diagn Cytopathol* 2011, 39: 313-317
- 20) Ordonez, N. G.: The diagnostic utility of immunohistochemistry in distinguishing between mesothelioma and renal cell carcinoma: a comparative study. *Hum Pathol* 2004, 35: 697-710
- 21) Butnor, K. J., Nicholson, A. G., Allred, D. C., et al.: Expression of renal cell carcinoma-associated markers erythropoietin, CD10, and renal cell carcinoma marker in diffuse malignant mesothelioma and metastatic renal cell carcinoma. *Arch Pathol Lab Med* 2006, 130: 823-827
- 22) Miettinen, M., McCue, P. A., Sarlomo-Rikala, M., et al.: GATA3: a multispecific but potentially useful marker in surgical pathology: a systematic analysis of 2500 epithelial and nonepithelial tumors. *Am J Surg Pathol* 2014, 38: 13-22
- 23) Matsuyama, A., Hisaoka, M., Iwasaki, M., et al.: TLE1 expression in malignant mesothelioma. *Virchows Arch* 2010, 457: 577-583
- 24) Chapel, D. B., Husain, A. N., Krausz, T., et al.: PAX8 Expression in a Subset of Malignant Peritoneal Mesotheliomas and Benign Mesothelium has Diagnostic Implications in the Differential Diagnosis of Ovarian Serous Carcinoma. *Am J Surg Pathol* 2017, 41: 1675-1682
- 25) Wu, D., Hiroshima, K., Matsumoto, S., et al.: Diagnostic usefulness of p16/CDKN2A FISH in distinguishing between sarcomatoid mesothelioma and fibrous pleuritis. *Am J Clin Pathol* 2013, 139: 39-46
- 26) Cigognetti, M., Lonardi, S., Fisogni, S., et al.: BAP1 (BRCA1-associated protein 1) is a highly specific marker for differentiating mesothelioma from reactive mesothelial proliferations. *Mod Pathol* 2015, 28: 1043-1057
- 27) Hida, T., Hamasaki, M., Matsumoto, S., et al.: Immunohistochemical detection of MTAP and BAP1 protein loss for mesothelioma diagnosis: Comparison with 9p21 FISH and BAP1 immunohistochemistry. *Lung Cancer* 2017, 104: 98-105
- 28) Wu, D., Hiroshima, K., Yusa, T., et al.: Usefulness of p16/CDKN2A fluorescence in situ hybridization and BAP1 immunohistochemistry for the diagnosis of biphasic mesothelioma. *Ann Diagn Pathol* 2017, 26: 31-37
- 29) Churg, A., Colby, T. V., Cagle, P., et al.: The separation of benign and malignant mesothelial proliferations. *Am J Surg Pathol* 2000, 24: 1183-1200
- 30) 環境省：平成28年度石綿関連疾患に係る医学的所見の解析調査業務（FISH法を用いた中皮腫診断法の開発に関する調査編）報告書（主任研究者：廣島健三）。（2016）。
- 31) Tochigi, N., Attanoos, R., Chirieac, L. R., et al.: p16 Deletion in sarcomatoid tumors of the lung and pleura. *Arch Pathol Lab Med* 2013, 137: 632-636

- 32) Hwang, H. C., Pyott, S., Rodriguez, S., et al.: BAP1 Immunohistochemistry and p16 FISH in the Diagnosis of Sarcomatous and Desmoplastic Mesotheliomas. *Am J Surg Pathol* 2016, 40: 714-718
- 33) Galateau Salle, F., Le Stang, N., Tirode, F., et al.: Comprehensive Molecular and Pathologic Evaluation of Transitional Mesothelioma Assisted by Deep Learning Approach: A Multi-Institutional Study of the International Mesothelioma Panel from the MESOPATH Reference Center. *J Thorac Oncol* 2020, 15: 1037-1053
- 34) Baranov, E., McBride, M. J., Bellizzi, A. M., et al.: A Novel SS18-SSX Fusion-specific Antibody for the Diagnosis of Synovial Sarcoma. *Am J Surg Pathol* 2020, 44: 922-933
- 35) Nakashima, Y., Inamura, K., Ninomiya, H., et al.: Frequent expression of conventional endothelial markers in pleural mesothelioma: usefulness of claudin-5 as well as combined traditional markers to distinguish mesothelioma from angiosarcoma. *Lung Cancer* 2020, 148: 20-27
- 36) Marchevsky, A. M., LeStang, N., Hiroshima, K., et al.: The differential diagnosis between pleural sarcomatoid mesothelioma and spindle cell/pleomorphic (sarcomatoid) carcinomas of the lung: evidence-based guidelines from the International Mesothelioma Panel and the MESOPATH National Reference Center. *Hum Pathol* 2017, 67: 160-168
- 37) Berg, K. B., & Churg, A.: GATA3 Immunohistochemistry for Distinguishing Sarcomatoid and Desmoplastic Mesothelioma From Sarcomatoid Carcinoma of the Lung. *Am J Surg Pathol* 2017, 41: 1221-1225
- 38) Amatya, V. J., Kushitani, K., Mawas, A. S., et al.: MUC4, a novel immunohistochemical marker identified by gene expression profiling, differentiates pleural sarcomatoid mesothelioma from lung sarcomatoid carcinoma. *Mod Pathol* 2017, 30: 672-681
- 39) Berg, K. B., Ionescu, D., & Churg, A. M.: MUC4 Staining in Sarcomatoid Carcinomas. *Mod Pathol* 2019, 32: 157
- 40) Terra, S., Roden, A., Yi, E., et al.: Utility of immunohistochemistry for MUC4 and GATA3 to aid in the distinction of sarcomatoid mesothelioma from sarcomatoid carcinoma. *Mod Pathol Supplement 2*. 2019, 32: 59
- 41) Whitaker, D., Henderson, D. W., & Shilkin, K. B.: The concept of mesothelioma in situ: implications for diagnosis and histogenesis. *Semin Diagn Pathol* 1992, 9: 151-161
- 42) Henderson, D. W., Shilkin, K. B., & Whitaker, D.: Reactive mesothelial hyperplasia vs mesothelioma, including mesothelioma in situ: a brief review. *Am J Clin Pathol* 1998, 110: 397-404
- 43) Churg, A., Hwang, H., Tan, L., et al.: Malignant mesothelioma in situ. *Histopathology* 2018, 72: 1033-1038
- 44) Churg, A., Galateau-Salle, F., Roden, A. C., et al.: Malignant mesothelioma in situ: morphologic features and clinical outcome. *Mod Pathol* 2020, 33: 297-302
- 45) Hasegawa, S., Kondo, N., Matsumoto, S., et al.: Practical approaches to diagnose and treat for T0 malignant pleural mesothelioma: a proposal for diagnostic total parietal pleurectomy. *Int J Clin Oncol* 2012, 17: 33-39
- 46) Stevers, M., Rabban, J. T., Garg, K., et al.: Well-differentiated papillary mesothelioma of the peritoneum is genetically defined by mutually exclusive mutations in TRAF7 and CDC42. *Mod Pathol* 2019, 32: 88-99
- 47) Offin, M., Aguirre, N., Yang, S. R., et al.: Clinical Characteristics and Outcomes of Patients with Well-Differentiated Papillary Peritoneal Mesothelial Tumors. *Ann Surg Oncol* 2024, 31: 7973-7977
- 48) Shrestha, R., Nabavi, N., Volik, S., et al.: Well-Differentiated Papillary Mesothelioma of the Peritoneum Is Genetically Distinct from Malignant Mesothelioma. *Cancers (Basel)* 2020, 12
- 49) Lee, H. E., Molina, J. R., Sukov, W. R., et al.: BAP1 loss is unusual in well-differentiated papillary mesothelioma and may predict development of malignant mesothelioma. *Hum Pathol* 2018, 79: 168-176
- 50) Hjerpe, A., Ascoli, V., Bedrossian, C. W., et al.: Guidelines for the cytopathologic diagnosis of epithelioid and mixed-type malignant mesothelioma. Complementary statement from the international mesothelioma interest group, also endorsed by the international academy of cytology and the papanicolaou society of cytopathology. *Acta Cytol* 2015, 59: 2-16
- 51) Hwang, H. C., Sheffield, B. S., Rodriguez, S., et al.: Utility of BAP1 immunohistochemistry and p16 (CDKN2A) FISH in the diagnosis of malignant mesothelioma in effusion cytology specimens. *Am J Surg Pathol* 2016, 40: 120-126
- 52) Walts, A. E., Hiroshima, K., McGregor, S. M., et al.: BAP1 immunostain and CDKN2A (p16) FISH analysis: clinical applicability for the diagnosis of malignant mesothelioma in effusions. *Diagn Cytopathol* 2016, 44: 599-606
- 53) Hiroshima, K., Wu, D., Hasegawa, M., et al.: Cytologic differential diagnosis of malignant mesothelioma and reactive mesothelial cells with FISH analysis of p16. *Diagn Cytopathol* 2016, 44: 591-598
- 54) Scherpereel, A., Astoul, P., Baas, P., et al.: Guidelines of the European Respiratory Society and the European Society of Thoracic Surgeons for the management of malignant pleural mesothelioma. *Eur Respir J* 2010, 35: 479-495
- 55) BTS statement on malignant mesothelioma in the UK, 2007. *Thorax* 2007, 62 Suppl 2: ii1-ii19.

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- 1) Husain AN et al., Guidelines for Pathologic Diagnosis of Mesothelioma: 2023 Update of the Consensus Statement From the International Mesothelioma Interest Group. Arch Pathol Lab Med. 2024;148(11):1251-1271.
- 2) Nabeshima K, Hamasaki M, Furukawa T, Matsumoto S, Takizawa K, Adachi M, Goto Y. Utility and limitations of Merlin immunohistochemistry for mesothelioma diagnosis in tissue sections and cell blocks. Ann Diagn Pathol. 2025 Aug;77:152475.
- 3) Hida T, Hamasaki M, Matsumoto S, Sato A, Tsujimura T, Kawahara K, Iwasaki A, Okamoto T, Oda Y, Honda H, Nabeshima K. BAP1 immunohistochemistry and p16 FISH results in combination provide higher confidence in malignant pleural mesothelioma diagnosis: ROC analysis of the two tests. Pathol Int. 2016 Oct;66(10):563-570.
- 4) Yoshimura M, Kinoshita Y, Hamasaki M, Matsumoto S, Hida T, Oda Y, Iwasaki A, Nabeshima K. Highly expressed EZH2 in combination with BAP1 and MTAP loss, as detected by immunohistochemistry, is useful for differentiating malignant pleural mesothelioma from reactive mesothelial hyperplasia. Lung Cancer. 2019 Apr;130:187-193.
- 5) Sa-Ngiamwibool P, Hamasaki M, Kinoshita Y, Matsumoto S, Sato A, Tsujimura T, Kasai T, Hiroshima K, Kushitani K, Takeshima Y, Kawahara K, Iwasaki A, Nabeshima K. Challenges and limitation of MTAP immunohistochemistry in diagnosing desmoplastic mesothelioma/sarcomatoid pleural mesothelioma with desmoplastic features. Ann Diagn Pathol. 2022 Oct;60:152004.

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